

**кольца для верхних моляров различных производителей**

размер, мм	Dentaurum Standart	Dentaurum Dentaform	GAC Snap-Fit	ORMCO MARK II	ORMCO Standard	ORMCO ULTIMA®	ORMCO WASHBON®	TP LABS	UNITEK CONTOURED	UNITEK GENERAL PURPOSE	Leone WEB	SIA Molar Bands
7.94.	1		1									
8.06.	2											
8.18.	3		3				1					
8.30.	4			2	2							
8.42.	5			3	3		2			30		
8.54.	6		4							30+		
8.66.	7			5	5	5						
8.78	8	1	5	6		7	5	2		31		
8.89	9	2		7	8		6	3	32	32		
8.99		3										
9.02.	10		6	9	9			5				
9.09		4										29.5.
9.14.	11			10	11	11	9	7		33		
9.20		5									1	30.5
9.25.	12		7	11			10	8	33+	33+		
9.30		6									2	31
9.38.	13			12	12	12		9	34	34		
9.40		7									3	
												31.5
9.50	14	8			13	13		10			4	
9.61	15	9		14	15		12	12				
											5	
9.71		10										32
9.74.	16		9			15	13	13	35	35		
9.81		11									6	
9.86.	17			16		16	14	14	35+	35+		
9.92		12									7	32.5
9.98.	18		11		17	17	15	15				
10.02		13									8	33
10.10.	19			18	18	18		17		36+	9	

10.12		14									10	33.5
10.23	20	15	14	19	19	20		19	36+			34
10.33.	21	16	15	20	20	21		21			11	34.5
											12	
10.43		17									13	35
10.46.	22			21	21							
10.54		18										
10.58.	23		17		22		20	23		38	14	
10.64		19									15	35.5
10.70.	24		18	22	23	23	21	24	38			
10.74		20									16	36
10.82.	25		19		24	24			38+	39		
10.85		21									17	
											18	36.5
10.95.	26	22		24	25		24	27				
											19	
11.05	27	23		26	27	27	27	30		39+	20	37
11.16		24									21	
11.18.	28		24			29	28		40	40+		
11.26		25										37.5
11.30.	29			28	29			32	40+		22	
11.36		26										
11.42.	30		26	29	30	31	29					
11.47		27										38
11.54.	31			30	31	32	30		41		23	
11.57		28										
11.67	32	29	28	32						41+	24	38.5

Эксклюзивный представитель компании «Дентаурум» в Украине

СП «Промед», ул. Лысенка 4а, Киев, 01034, Украина

Тел: (044) 278-7339, 278-7346, 531-9123, (050) 441-8554

Факс: (044) 278-6245

E-mail: [admin@promed.ua](mailto:admin@promed.ua) [www.promed.ua](http://www.promed.ua)